

Patient Information	Specimen Information	Client Information

COMMENTS:

Test Name	In Range	Out Of Range	Reference Range	Lab
PARTIAL THROMBOPLASTIN TIME, ACTIVATED	26		23-32 sec	

This test has not been validated for monitoring unfractionated heparin therapy. For testing that is validated for this type of therapy, please refer to the Heparin Anti-Xa assay (test code 30292).

For additional information, please refer to <http://education.QuestDiagnostics.com/faq/FAQ159> (This link is being provided for informational/educational purposes only.)

PERFORMING SITE:

SPECIMEN: